

# PREempt™ CS20 Tool & Implement Sterilant

Sterilizes in 20 minutes with no dilution required  
and reusable for up to 14 days.

## POWERFUL

Kills spores and  
fungi on tools in  
20 MINUTES.

## APPROVED

Health Canada registered  
high level disinfectant.  
Class II Medical Device.

## COMPATIBLE

With your professional  
stainless steel and  
hard plastic tools and  
implements.



## CONVENIENT

No diluting required.  
Re-use up to 14 days  
in a tool soaker.<sup>1</sup>

## GENTLE

No fumes or special  
ventilation required.

## RESPONSIBLE

Biodegradable. Active  
ingredient breaks down  
into water and oxygen.

<sup>1</sup> High-level disinfection is recommended for semi-critical tools or implements, that come in contact with mucous membranes or non-intact skin during use, but do not ordinarily penetrate the skin or enter normally sterile areas of the body. In personal service settings, this can include tools such as reusable foot files, tweezers used for ingrown hairs, microblade handles, nail clippers and cuticle pushers, or needle holders. Health Canada is responsible for oversight of high-level disinfectant label claims and has approved the re-use of PREempt™ CS20 for 14 days. Exceptions to this allowance occur if the solution is visibly soiled, or if the indicator strip fails to produce the required measurement. In these cases, the solution should be discarded even if the 14-day period has not yet been reached.